

## Kiln

466965

### Maintenance Code

W/B+Steam - Clean with water-based cleanser, diluted household bleach, or steam clean.



### General Recommendations

Treat stains promptly. The longer a stain remains on the textile, the more difficult it will be to remove. A professional furniture cleaning service is recommended. Vacuum regularly using the proper attachment to avoid pilling and/or clean with water-based cleanser, diluted household bleach, or steam as needed. Always pre-test any cleaning method in an inconspicuous area. Several light applications of a cleaning agent are less harmful than a single concentrated application. Repeat stain treatment several times if necessary.

Do not use solvent-based cleaning agents such as lacquer thinners, nail polish remover, or pine oil. These agents may cause discoloration over time.

Do not over saturate the textile or let water pool in seams.

### Stain Treatment - Water-Based Cleansers

- For oil and water-based stains, use a clean, soft, white cloth, a natural sponge, or a nylon soft-bristle brush with lukewarm, soapy water. Brush the stain with light, quick strokes. Avoid a scrubbing motion and work from the edge of the stain toward the center.
- For tougher oil-based stains, dab a liberal amount of rubbing alcohol on the stain and let it sit for one minute before applying a water based cleanser. Confirm the construction of your textile is bleach cleanable. If so, and additional stain treatment is needed, use a 10% solution of diluted household bleach. Thoroughly rinse with cool, clean water and vacuum extract the textile to remove any residue.
- For solid stains (e.g. mud) allow the stain to dry thoroughly and gently break up any hardened areas. Use a nylon, soft bristle brush to remove all loose particles.
- For insoluble stains (e.g. chewing gum, candle wax, etc.), use a plastic bag filled with ice to rub the stain until the substance becomes brittle. Use a clean, dull, flat knife to gently remove the substance. If the stain persists, blot the area with a clean, soft, white cloth, or natural sponge and lukewarm soapy water.
- For blood and bodily fluid stains, follow the oil-based stain treatment above. If the stain persists, gently dab a pH-neutral cleanser, upholstery spot cleaner, or an equal mixture of water and enzyme detergent (e.g. laundry detergent) on the stain and let it sit for one minute.

Rinse all stain treatments thoroughly with a natural sponge moistened slightly with cool, clean water and completely remove all cleaning solution residue as left over residue can attract soil over time. Extra care should be taken when cleaning dark colors of woven textiles to avoid creating light marks that can be caused by color transfer. Avoid excess rubbing to prevent discoloration and damage to the face of the material.

Allow the textile to dry thoroughly and vacuum before using.

### Stain Treatment - Steam Cleaning

Before cleaning, fill steam cleaner with fresh, clear water. Pre-treat stains with a professional steam cleaning pre-treatment solution or an equal mixture of enzyme powder detergent and water.

Kiln

466965

## Healthcare Recommendations

Healthcare facilities should follow safety protocols for biological hazards and/or anti-viral treatments during the cleaning process. For a complete list of approved hospital grade cleansers, contact your sales representative.

## Note

Maharam can neither guarantee, nor take responsibility for specific results. Under no circumstances shall Maharam be liable for any incidental, consequential, or other damage resulting from negligence, misuse, mishandling, or exposure of any chemicals used in connection with maintenance. Information relating to chemicals/cleaning solvents, environmental use, and human exposure should be reviewed with appropriate suppliers.

# Maintenance Guidelines

## Kiln

466965

### Approved Cleaners

#### 3M

3M HB Quat Disinfectant Cleaner Concentrate  
3M HB Quat Disinfectant Cleaner RTU  
3M Neutral Quat Disinfectant Cleaner Concentrat  
3M Quat Disinfectant Cleaner Concentrate

#### Amsan

Sureblend Hospital Disinfectant Cleaner

#### Brighton Professional

Enzyme Plus Odor Eliminator Deodorizer

#### Buckeye International

Sanicare Quat 256

#### Core Products

Unbelievable! Stain Remover and Deodorizer

#### Diversey

General Purpose Cleaner #15  
Glance NA Glass Multi-Surface Cleaner Non-Amm  
Stride Floral Neutral Cleaner  
Stride Fragrance Free SC Neutral Cleaner  
Virex II 256

#### Ecolab

A-456-N  
Eco-Star Detergent  
Enverros Fiberfresh Ultra Pro  
Qc 31 Neutral Cleaner  
Viva

#### Medline Industries

Micro-Kill Bleach Cermicidal Bleach Wipes

#### PDI

Sani-Cloth Bleach Germicidal Disposable Wipe

#### Professional Products

Quick Fill 120

#### Spartan Chemical Company

Super HDQ Neutral

#### Sustainable Earth by Staples

63 Concentrated Odor Eliminator  
64 Neutral Multi-Use Cleaner  
67 Carpet Spotter RTU

#### Sysco

Carpet Textile Stain Digester

#### The Clorox Company

Clorox Healthcare Bleach Germicidal Wipes

#### Various

1:10 Solution of Diluted Household Bleach  
Water-Based Cleaners

#### Waxie Enterprises

Solsta 700 Disinfectant Cleaner

#### Wood Wyant

Ed - Everyday Disinfectant  
Ultraquat Germicidal Detergent